

## **NEXT GENERATION SEQUENCING ANALYSIS(ASSIGNMENT PART2 )**

### **Data Retrieval (NCBI SRA) & Introduction to data types**

- Read Quality Check (FastQC & Cutadapt)
- Alignment of reads using reference Genome (BWA)
- Visualization of mapped reads (IGV)

### ***Paired End-RAW DATA***

Paired-end sequencing enables the sequencing of both ends of DNA fragments, producing high-quality, alignable sequence data. In this process, two files are generated: one containing all the forward sequences and the other containing the reverse sequences, with both files aligned in the same order. This method is particularly useful for detecting genomic rearrangements, repetitive sequence elements, gene fusions, and novel transcripts.

### **Treated – R1norm&R2norm**

### **Control – R1TUMOR&R2TUMOR**

The FASTQ format is a text-based method for storing biological sequences along with their corresponding quality scores. Each sequence in a FASTQ file spans four lines. The last line contains ASCII characters representing the quality scores of the bases in the second line. These quality scores, which identify nucleotide bases, are measured using the Phred quality score, indicating the probability of an error in base calling. The Phred scale's primary function is to automatically determine accurate, quality-based consensus sequences. For brevity, both the sequence letter and the quality score are encoded with ASCII characters.

### **FASTQC**

FastQC offers an easy method for performing quality control checks on raw sequence data from high-throughput sequencing pipelines. It is compatible with data produced by Illumina platforms.

Installation of FASTQC in LINUX

**sudo apt-get install fastqc**

To get the Quality report of raw data

**fastqc R1norm R2norm R1TUMOR R2TUMOR**

If the read quality is satisfactory, you can proceed with the mapping process. If not, use Cutadapt software to clean the data.

### **CUTADAPT**

Cutadapt identifies and eliminates adapter sequences, primers, poly-A tails, and other unwanted sequences from high-throughput sequencing reads.

To download and install Cutadapt

**sudo pip install cutadapt**

To trim a 3' adapter, with AACCGGTT as the adapter to be removed, reads are processed from input.fastq and the results are saved in output.fastq.

**cutadapt -a AACCGGTT -o cutout\_1.fastq read1.fastq**

To trim a 5' adapter, with AACCGGTT as the adapter to be removed, reads are processed from input.fastq and the results are saved in output.fastq.

**cutadapt -g AACCGGTT -o cutout\_1.fastq read1.fastq**

FOR paired end(R1norm&R2norm)

**cutadapt -q 10 -a ADAPTER --minimum-length 20 -o tmp.1.fastq -p tmp.2.fastq R1norm.fastq & R2norm.fastq**

**cutadapt -q 15 -a ADAPTER --minimum-length 20 -o trimmed.2.fastq -p trimmed.1.fastq tmp.2.fastq tmp.1.fastq**

FOR paired end(R1TUMOR&R2TUMOR)

**cutadapt -q 10 -a ADAPTER --minimum-length 20 -o tmp.1.fastq -p tmp.2.fastq R1TUMOR.fastq & R2TUMOR.fastq**

**cutadapt -q 15 -a ADAPTER --minimum-length 20 -o trimmedtumor.2.fastq -p trimmedtumor.1.fastq tmptumor.2.fastq tmptumor.1.fastq**

Again check the quality of reads; if everything is good, go ahead with mapping.

### **Mapping on Reference Genome using BWA.**

DOWNLOAD THE REFERENCE GENOME FROM UCSC GENOME BROWSER. T2T-CHM13:

First, prepare the 'hs1.fa.gz' file.

To analyze sequence data, aligning it to a reference genome is essential. BWA is a software package designed for mapping low-divergence sequences against large reference genomes, like the human genome.

BWA first needs to construct index for the reference genome.

```
gunzip hs1.fa.gz
bwa index -a bwtsv hs1.fa
```

```
sudo falloate -l 16G /swapfile
sudo chmod 600 /swapfile
sudo mkswap /swapfile
sudo swapon /swapfile
free -m
```


## Using Mem

To align the reads on the reference genome  
For paired end:

```
bwa mem hs1.fa trimmed.1.fastq trimmed.2.fastq > bwa_1.sam
bwa mem hs1.fa trimmedtumor.1.fastq trimmedtumor.2.fastq > bwa_2.sam
```

The **GALAXY** server was utilized to run Strelka or VarScan for the purpose of identifying somatic mutations.

The tumour and normal SAM files, along with the reference genome file, were uploaded to the GALAXY server.



With four output files in VCF format, we can perform various analysis.

```
Galaxy9-[Strelka_Somatic_on_data_1_data_4_and_data_3_SNVs_vcf]
Galaxy8-[Strelka_Somatic_on_data_1_data_4_and_data_3_Indels_vcf]
Galaxy6-[VarScan_somatic_indel_calls_on_data_4_data_3_and_data_1]
Galaxy5-[VarScan_somatic_SNP_calls_on_data_4_data_3_and_data_1]
```

**IGV (Integrative Genomics Viewer):** A powerful visualization tool that allows you to view VCF files and explore genetic variations in a genome browser interface given below.

